

Petri Dish Selection Guide

Manufactured in high capacity automated production lines resulting in high quality, value priced products!

100mm x 15mm Petri Dish Options



STACKABLE

Lids contain stacking rings for increased stability

229695



STACKABLE/SLIDEABLE

Features partial edge ring on base which allows a single dish to be slid from a large stack and enhances stability when stacked

229694




SLIPPABLE

Designed for automated filling and streaking systems, features side arrows and ISO Mark Target for automation

229697


Maximize time and resources with smaller size dishes and divided 100mm dishes:

Smaller Size Petri Dishes



229638, 229665

Compartment Petri Dishes




Molded dividers create separate compartments while allowing airflow

229682, 229683, 229684

Part No.	Product Description	Working Volume	Growth Area	Pack Description	Qty/Case	Dish Style
229638	35mm x 10mm Petri Dish	5mL	8.5cm ²	10/Bag	960	Stackable
229665	60mm x 15mm Petri Dish, Stackable	7-8mL	21cm ²	20/Bag	500	Stackable
229694	100mm x 15mm Petri Dish, Slideable	15-16mL	55cm ²	20/Bag	500	Stackable/Slideable
229695	100mm x 15mm Petri Dish, Stackable	15-16mL	55cm ²	20/Bag	500	Stackable
229697	100mm x 15mm Petri Dish, Slippable	15-16mL	55cm ²	20/Bag	500	Slippable
229682	100mm x 15mm Petri Dish, 2 Compartments	7-8mL	28.3cm ²	20/Bag	500	Compartment
229683	100mm x 15mm Petri Dish, 3 Compartments	5mL	18.9cm ²	20/Bag	500	Compartment
229684	100mm x 15mm Petri Dish, 4 Compartments	4mL	14.1cm ²	20/Bag	500	Compartment



Loops, Needles & Spreaders, oh my!

Popular accessories to make working with Petri Dishes simple and easy!



- Disposable, sterile tools eliminate the need to sterilize instruments between cultures
- L-Shaped Cell Spreader is ideal for evenly spreading bacterial cultures; angled tip prevents damaging agar surfaces
- Loops and Needles are flexible for easy access to small-volume tubes
- All tools are color-coded for easy identification
- Loops and L-Shaped Cell Spreader available individually wrapped in peel-to-open paper/plastic wrappers

L-SHAPED SPREADERS					AF, S, R/D F, NP, USP, PP
Part No.	Product Description	Color	Length	Pack Description	Qty/Case
229616	L-Shaped Spreader	Green	145mm	100 Individually Wrapped/Bag	100
229617	L-Shaped Spreader	Green	145mm	10/Pack	500
INOCULATING LOOPS AND NEEDLES					AF, S, R/D F, NP, USP, PP
Part No.	Product Description	Color	Length	Pack Description	Qty/Case
229601	1µL Inoculating Loop	White	218mm	20/Pack	2000
229602	1µL Inoculating Loop	White	218mm	100 Individually Wrapped/Bag	600
229610	10µL Inoculating Loop	Blue	221mm	20/Pack	2000
229612	10µL Inoculating Loop	Blue	221mm	100 Individually Wrapped/Bag	600
229615	Inoculating Needle	Yellow	221mm	20/Pack	2000



Lit Code: 45A

