


Check Out What's **NEW!**



Experience Easier Pipetting



Low-force Tip Solution Pipette Tips (LFTS)

- Soft plastic on mounting end allows for low-force tip acquisition
- Compatible with Low-force pipetting systems
- More sizes available

Part No.	Description	Pack	Qty/Case
229081	300µL LFTS Filter Pipette Tips	96/Rack	960
229082	1200µL LFTS Filter Pipette Tips	96/Rack	960
229083	300µL LFTS Pipette Tips	96/Rack	960
229084	1200µL LFTS Pipette Tips	96/Rack	960



The CELLTREAT® Difference

FEATURED PRODUCTS

Discover how square can save!

5mL Reagent Reservoirs

Don't risk losing expensive samples and solutions!



- Unique square, pyramid bottom reservoir enables easy and complete transfer of precious samples
- Eliminate contamination risk of repeat pipetting from a narrow tube
- Both miniature 5mL designs consume less plastic, supporting sustainability goals

Part No.	Description	Pack	Qty/Case
229288	5mL Reservoir, Individually Wrapped	1/Bag	80
229289	5mL Reservoir, Bulk	5/Bag	200
229292	5mL Reservoir, Square Pyramid Bottom, Ind Wrap	1/Bag	80
229293	5mL Reservoir, Square Pyramid Bottom, Bulk	5/Bag	200

The better bottle is here!

PETG Square Media Bottles

When it comes to maximizing limited refrigerated space, round is out and square is in!

- Safely store solutions in PETG Square Media Bottles - no risk of breaking glass
- PETG material reduces gas permeability to improve solution stability
- Sterile bottles available individually wrapped or bulk tray packed



Part No.	Description	Pack	Qty/Case	Part No.	Description	Pack	Qty/Case
230750	30mL Media Bottle, Square	1/Bag	24	230760	30mL Media Bottle, Square	24/Tray	48
230751	60mL Media Bottle, Square	1/Bag	24	230761	60mL Media Bottle, Square	24/Tray	48
230752	125mL Media Bottle, Square	1/Bag	24	230762	125mL Media Bottle, Square	12/Tray	24
230753	250mL Media Bottle, Square	1/Bag	24	230763	250mL Media Bottle, Square	12/Tray	24
230754	500mL Media Bottle, Square	1/Bag	24	230764	500mL Media Bottle, Square	12/Tray	24
230755	1000mL Media Bottle, Square	1/Bag	24	230765	1000mL Media Bottle, Square	12/Tray	24

Meet the Newest Members of the CELLTREAT Scientific Products & VistaLab Technologies Family

Protein processing that won't break the bank!

Centrifugal Filter Units (MWCO Tubes)

Swap your cap color; green is the new pink!

- Efficiently separate protein samples by molecular weight size exclusion (5, 10, 30, 50, and 100kDa)
- Cost savings: small case quantities are budget-friendly
- Protein purification, buffer exchange, sample concentration, and more!
- Dual vertical PES filters efficiently process samples with >90% sample recovery



Part No.	Description	Pack	Qty/Case
230708	4mL Centrifugal Filter, MWCO 5kDa	8/Box	8
230710	4mL Centrifugal Filter, MWCO 10kDa	8/Box	8
230712	4mL Centrifugal Filter, MWCO 30kDa	8/Box	8
230714	4mL Centrifugal Filter, MWCO 50kDa	8/Box	8
230716	4mL Centrifugal Filter, MWCO 100kDa	8/Box	8

Part No.	Description	Pack	Qty/Case
230720	15mL Centrifugal Filter, MWCO 5kDa	8/Box	8
230722	15mL Centrifugal Filter, MWCO 10kDa	8/Box	8
230724	15mL Centrifugal Filter, MWCO 30kDa	8/Box	8
230726	15mL Centrifugal Filter, MWCO 50kDa	8/Box	8
230728	15mL Centrifugal Filter, MWCO 100kDa	8/Box	8

Add power to your research!

Permeable Cell Culture Inserts

NEW! 0.1µm Membranes

- Easily mimic *in vivo* conditions
- Assay versatility: cell-cell interactions, drug transport, angiogenesis, migration, invasion, and more!
- 0.4µm, 3.0µm, and 8.0µm systems also available, see all products below



	6 Well		12 Well		24 Well	
	PC	PET	PC	PET	PC	PET
0.1µm	230600	230606	230614	230620	230628	230634
0.4µm	230601	230607	230615	230621	230629	230635
3.0µm	230603	230609	230617	230623	230631	230637
8.0µm	230605	230611	230619	230625	230633	230639

Low retention, high performance!

Universal Fit Pipette Tips

NEW! 10µL Extended Length Reload System

NEW! 1250µL Racked Filter Tips

- Reduce plastic waste and save time and money by refilling empty racks
- 1250µL racked filter tips are the perfect companion for large volume multichannel pipettes
- Full line of filter and non-filter pipette tips available at www.celltreat.com



Part No.	Description – Reload System	Pack	Qty/Case
229053	10µL Extended Length Reload System	960 Tips/Reload	1920

Part No.	Description – Racked Filter Tips	Pack	Qty/Case
229023	1250µL Filter Pipette Tips	96/Rack	960


Repeat Pipetting is Faster and Easier with ali-Q™ 2



Compatible with ANY serological pipet!



Set your ideal speed for any solution or application



Turn the dial to set your exact aliquot volume



Benefits

- ✓ **SAVE TIME**
No need to watch the pipet graduation lines
- ✓ **SAVE MONEY**
Fewer repeated experiments due to pipetting errors
- ✓ **PEACE OF MIND**
Take the guesswork out of repeat pipetting!



LEARN MORE!
Watch the video and discover the difference.



ALIQOT:
equal repeat dispensing!



Press and hold aliquot button for accurate dispensing

Model	Part No.	Aliquoting Speed	Aliquoting Range
ali-Q™ 2 VS, 5mL	2100-1005	Variable	0.5-5mL
ali-Q™ 2 VS, 10mL NEW!	2100-1010	Variable	1.0-10mL

Trade-In Trade-Up!

YOU DESERVE AN UPGRADE!

Trade in your current pipet controller and save on ali-Q 2!

☎ 877.345.4620 ✉ service@celltreat.com
 🌐 www.celltreat.com/pipette-trade-in



Lit Code: 42A



CELLTREAT
SCIENTIFIC PRODUCTS
www.celltreat.com



VistaLab
TECHNOLOGIES
A subsidiary of CELLTREAT Scientific Products