

PIPETTE TIP ADAPTER COMPATIBILITY GUIDE

PIPETTE BRAND/ MANUFACTURER		EPPENDORF RESEARCH (S)					EPPENDORF RESEARCH (MULTI)					EPPENDORF RESEARCH PLUS			GILSON PIPETMAN (SINGLE)					MULTI	THERMO FINNPIPETTES (SINGLE)				SARTORIUS BIOHIT PROLINE			RAININ PIPET LITE XLS												
MODEL/RANGE		0.1-2.5µL	0.5-10µL	2-20µL	10-100µL	20-200µL	100-1000µL	0.5-105µL	5-100µL	50-100µL	20-300µL	30-300µL	50-1200µL	0.5-2.5µL	5-100µL	2-200µL	50-1000µL	P2, 0.1-2µL	P10, 0.2-10µL	P20, 2-20µL	P100, 10-100µL	P200, 20-200µL	P1000, 100-1000µL	P200 M8, 20-200µL	40-200µL	100-1000µL	200-1000µL	5-50µL	50-300µL	0.1-2.5µL	0.5-10µL	20-200µL	100-1000µL	SL-2, 0.1-2µL	SL-10, 0.5-10µL	20-200µL	100-1000µL			
VOLUME	PART NO.	UNIVERSAL FIT PIPETTE TIPS																																						
10µL	229042	✓					✓											✓	✓																✓	✓				
	229052		✓					✓											✓	✓																✓	✓			
	229032		✓						✓											✓	✓															✓	✓			
10µL, Extended	229043	✓	✓					✓						✓					✓																		✓			
	229033	✓	✓						✓					✓						✓																	✓			
200µL	229044				✓	✓			✓	✓	✓	✓			✓	✓				✓	✓	✓	✓	✓														✓		
	229054				✓	✓			✓	✓	✓	✓			✓	✓					✓	✓	✓	✓	✓													✓		
	229034				✓	✓			✓	✓	✓	✓			✓	✓					✓	✓	✓	✓	✓													✓		
300µL	229055				✓	✓			✓	✓	✓	✓			✓	✓					✓	✓	✓	✓	✓													✓		
	229035				✓	✓			✓	✓	✓	✓			✓	✓					✓	✓	✓	✓	✓													✓		
1000µL	229056						✓					✓											✓			✓	✓										✓		✓	
	229037							✓				✓											✓			✓	✓										✓		✓	
	229036								✓				✓										✓			✓	✓											✓		✓
VOLUME	PART NO.	UNIVERSAL FIT PIPETTE TIPS - FILTER																																						
10µL	229015	✓	✓					✓											✓	✓																✓	✓			
10µL, Extended	229016		✓						✓											✓																		✓		
20µL	229017				✓	✓			✓	✓	✓	✓			✓	✓					✓	✓	✓	✓	✓													✓		
100µL	229018				✓	✓			✓	✓	✓	✓			✓	✓					✓	✓	✓	✓	✓													✓		
200µL	229019				✓	✓			✓	✓	✓	✓			✓	✓					✓	✓	✓	✓	✓													✓		
300µL	229020				✓	✓			✓	✓	✓	✓			✓	✓					✓	✓	✓	✓	✓													✓		
1000µL	229022						✓					✓											✓			✓	✓											✓		✓
	229021							✓				✓											✓			✓	✓											✓		✓



LOOKING FOR VARIETY?

CELLTREAT pipette tips are available in a variety of sizes and options - bulk packed, low retention, low retention with filter, and in a reload system

NEED HELP ORDERING?

For order assistance, please email service@celltreat.com

WANT TO TRY A SAMPLE?

To request samples, please email samples@celltreat.com or complete the online form at: bit.ly/requestsamplenow

